

Instructions for Use

Willow Sync™
In-Bra Breast Pump

willowsync™

Instructions for Use are updated periodically, so
please check onewillow.com/SyncIFU for current version.



TABLE OF CONTENTS

Overview	02
Warnings	
Precautions	
Symbol definitions	
Indications for use	
Usage & operating life	
Product description	
Product overview	
Product parts	
Pump controls	
Usage Instructions	11
1. Getting started: Disassembly	
2. Getting started: Before first use	
3. Assembly	
4. Putting on the pump	
5. Pumping	
6. Transferring milk	
7. Cleaning	
8. Charging	
9. Troubleshooting	
Technical Information	35

OVERVIEW

WARNINGS

1. Do not operate where aerosol spray products are used or oxygen is being administered.
2. Use Product only for its intended use as described herein.
3. This Product is for a single user and should not be shared between users.
4. Only use the Product parts. Do not substitute or modify pump parts.
5. The Product contains parts, such as the USB Charging Cable, which may present a hazard. Keep out of reach of children.
6. Do not attempt to open or repair Motor or Receiver. No modification of this equipment is allowed. The Motor contains an internal non-replaceable lithium-ion battery pack.
7. The Motor is not a serviceable part.
8. Never operate the Product if USB Charging Cable is damaged or not working properly.
9. To protect against fire, electric shock and injury to users, do not immerse Motor, USB Cable, or Charger in water. If these components are immersed in water, do not operate the Product.
10. Do not operate if any component in the Product is damaged or not working properly.
11. Misuse of the Product may result in the risk of fire, electric shock and/or personal injury.
12. The Product is not provided sterile and must be cleaned, and all washable parts must be boiled for 5 minutes before first use.

-
13. If the USB Charging Cable casing or wiring becomes loose, separated, or frayed, stop use of the USB Charging Cable immediately and contact Willow. Stop using the Product immediately if you notice any smoke or burning from the Product.
 14. Do not use the Product while driving, when sleep/drowsy, or while lying down.
 15. Care should be taken when storing the Product. Improper storage (e.g., storing wet components, inappropriate temperatures, exposure to sunlight, etc.) may result in damage and cause harm to the user.
 16. The Product is not intended to be used during air transport.
 17. If you are a mother who is infected with Hepatitis B, Hepatitis C or Human Immunodeficiency Virus (HIV), pumping breast milk will not reduce or remove the risk of transmitting the virus to your baby through your breast milk.
 18. Monitor milk level on Container Bowls during pumping to ensure overflow does not occur.
 19. If you feel excessive warmth from the Product while pumping, stop pumping immediately and remove pump. Contact Customer Care.
 20. If irritation or discomfort occurs, discontinue use and see a doctor.
 21. Splashguards provided in this pump kit are only compatible with Willow Sync™ Containers.
 22. The Willow Sync Breast Pump is only compatible with Willow Sync Containers. Do not use Willow Go containers with the Willow Sync Breast Pump.

PRECAUTIONS

1. If package appears damaged or tampered with before use, do not use and contact Customer Care at 1-888-WILLOW-3.
2. If pumping is uncomfortable, change suction, turn off and reposition, or remove and discontinue session. Consult your healthcare professional.
3. Before each use visually inspect each part for signs of wear and tear (cracks, chips, tears, discoloration or deterioration). In the event that damage to the Product is observed, discontinue use until the parts have been replaced.
4. If discomfort is felt, then suction can be broken by inserting a finger between the breast and the flange.
5. Before use, ensure that your hands and the system are clean.
6. Read all instructions before using the Product.

SYMBOL DEFINITIONS

For definition of symbols, please refer to the Symbols Glossary on pages 46–47.

INDICATIONS FOR USE

The Product is a powered breast pump to be used by lactating women to express and collect milk from their breasts. The Product is intended for a single user.

CONTRAINDICATION - No known contraindications.

USAGE & OPERATING LIFE


The Willow Sync™ Breast Pump has an expected operating life of 12 months. The Pump Diaphragm, Container Assembly, and Flange are reusable parts and should be cleaned after each use. The expected operating life of these components is 3 months.

PRODUCT DESCRIPTION

The Willow Sync™ Breast Pump is a personal use, wearable breast pump for single and double pumping. The Breast Pump is not provided sterile. Prior to first use, wipe down Motor, Receiver, and Tubing with disinfectant wipes, such as 70-90% ethyl or isopropyl alcohol. Boil all washable parts (Pump Diaphragm, Duckbill Valves, Flanges, Sizing Inserts, Container Backs, Container Back Gaskets, Container Bowls, and Splashguards) for 5 minutes. For routine cleaning, see instructions provided on pages 29–32.

For replacement parts visit onewillow.com or contact Customer Care at 1-888-WILLOW-3

PRODUCT OVERVIEW



Pump


Controls everything for your
pumping session

Pump Cup


Part Name	Quantity
Motor	1
Receiver	1
Tubing	1
Flange (24mm)	2
Sizing Insert (21mm)	2
Pump Diaphragm	2
Duckbill Valve	2
Container Back	2
Container Back Gasket	2
Container Bowl, 5oz	2
Splashguard	2
USB Charging Cable	1

NON-WASHABLE PARTS

⚠ Do not submerge or run under water



Motor



Receiver



Tubing



USB Charging Cable

*USB Wall Charger
sold separately




NON-WASHABLE PARTS


WASHABLE PARTS



Sizing Insert
21 mm



Flange
24 mm




Pump Diaphragm



WASHABLE PARTS

WASHABLE PARTS



Container Back



Container Back
Gasket



Duckbill Valve



Container Bowl
5 ounces (150 mL)



Splashguard*


*The Splashguard does not serve as a leakproof cap. The container has a necessary venthole so air can escape while pumping.



CONTAINER WASHABLE PARTS


CONTROLS

The controls on the Pump will sync to the Pump Cup automatically




USAGE INSTRUCTIONS


1. GETTING STARTED: DISASSEMBLY



Remove Flange



Unhook Container



Remove Diaphragm

1. GETTING STARTED: DISASSEMBLY

4



Remove Container Bowl
from Container Back

5




Open Splashguard and un-loop around spout

*The Splashguard does not serve as a leakproof cap.
The container has a necessary vent hole so air can
escape while pumping.


1. GETTING STARTED: DISASSEMBLY

6



Remove Duckbill Valve
from Container Back

7



Remove Gasket from
Container Back

2. GETTING STARTED: BEFORE FIRST USE


Clean and Sanitize Parts

The Willow Sync™ Breast Pump is not provided sterile. Prior to first use, clean the washable and non-washable parts. For routine cleaning, see instructions provided on pages 29–32.



WASHABLE PARTS

Boil all washable parts for 5 minutes.



CONTAINER WASHABLE PARTS


2. GETTING STARTED: BEFORE FIRST USE




NON-WASHABLE PARTS

Wipe down Motor, Receiver, and Tubing with disinfectant wipes (such as 70-90% ethyl or isopropyl alcohol).

⚠ **Do not submerge or run under water**



Motor



Receiver



Tubing



NON-WASHABLE PARTS


2. GETTING STARTED: BEFORE FIRST USE

Check your Flange fit

Use the Nipple Sizing Tool included in the box to measure your nipple and find your Flange size. Proper fit leads to greater comfort. 21 mm and 24 mm included. Additional 13 mm, 15 mm, 17 mm, 19 mm and 27 mm sizes are available for purchase.


Plug in Motor until indicator light turns on

The Motor comes partially charged, but needs to be woken up before first use. Plug the USB charging cable into the Motor, then into USB wall charger (not included). Once woken up, the LED Indicator will pulse.




3. ASSEMBLY


Note: Ensure gasket is assembled correctly - improper assembly can cause leaks.



3. ASSEMBLY




Press Duckbill Valve onto
Container Back, aligning tab



Close Container Bowl
onto Container Back

3. ASSEMBLY

4



Loop Splashguard around spout, fold over, and press into spout.


*The Splashguard does not serve as a leakproof cap.
The container has a necessary vent hole so air can escape while pumping.

5




Place Diaphragm into Motor by aligning tabs.
Pull Diaphragm through the hole, so that it lies flat.

3. ASSEMBLY



Hook and press assembled Container onto Pump



Press Flange into assembled Pump

3. ASSEMBLY


For Double-Sided Pumping

8



Repeat steps 1-7 for Pump Cup

9



Plug Tubing into Pump and Pump Cup
(a slight gap is OK)

For Single-Sided Pumping

10




Close Plug into Tubing Port on Motor

4. ALIGNING THE PUMP

Bring Pump or Pump Cup to your breast, centering your nipple in the Flange tunnel.

Important: Correct nipple alignment is key to comfort, output and performance.


- Nipple must be positioned correctly in Flange tunnel before pumping to avoid discomfort.
- If you still experience discomfort with proper placement, you may need a different sizing insert. Use the [Nipple Sizing Tool](#) to help select the appropriate Flange size.
- Use a full-coverage nursing bra with nursing flap that supports the Pump or Pump Cup firmly at top and bottom. For best support, you do not need to go up a cup size.



Align nipple in the center of the Flange tunnel



Bring Pump or Pump Cup straight onto breast




Close nursing bra. While pumping, keep your upper body in an upright position, whether sitting or standing.

4. ALIGNING THE PUMP


Double Pumping

The Pump can be worn on either side.



Single Pumping

Keep Tubing Port closed.



Note: If you prefer the Tubing shorter, remove the Tubing Connectors and cut Tubing to desired length. Re-install Tubing Connectors by pushing them into each end of the Tubing.


5. PUMPING

Hold the Power button on Motor for one second to start and stop pumping. The LED Indicator will display green, orange or red for battery status, then turn off after 5 seconds for discreet pumping.

Pumping begins in Stimulation Mode with quick suckling. When you feel let-down, you may press the MODE button to switch to slower and stronger Expression Mode. The Pump automatically transitions to Expression Mode after 2 minutes.

Press the MODE button at any time to change the speed from slow to fast. Adjust suction strength with the Increase / Decrease buttons to find the combination that is most comfortable.

Note: Motor will not work when connected to a charger. Always disconnect Motor from charger before using. LED Indicator will be off during pumping unless a notification is triggered.



Color of light indicates battery status:


When Power button is pressed	While charging	
 Green Solid	 Green Pulsing	Enough for at least 3 pumping sessions (approximately 25 minutes) LED Indicator will turn off once battery is fully charged.
 Orange Solid	 Orange Pulsing	Enough for at least 1 pumping session (approximately 25 minutes)
 Red Solid	 Red Pulsing	Low battery - inadequately charged for session

	Stimulation Mode	Expression Mode
Number of Levels	10	10
Starting Level	1	Level used in Stimulation

Important: Avoid setting suction at a level that causes discomfort. If you feel excessive discomfort or pain, stop pumping, remove Pump and realign or discontinue session.

6. ENDING A SESSION AND TRANSFERRING MILK

1. When finished pumping, power off the Pump and disconnect the tubing, if using. Remove the Pump and Pump Cup from your nursing bra.
2. Transfer the milk to a storage bottle / bag before removing the Container from the Pump or Pump Cup. Do not wait before transferring the milk out of the container to avoid any backflow from the flange.



BREAST MILK STORAGE GUIDELINES

	Storage Location and Temperatures			
	Room Temperature 60–77° F	Refrigerator Below 39.2° F	Freezer Below 24.8° F	Insulated Cooler Below 59° F
Freshly expressed or pumped	Up to 4 hours optimal 6–8 hours acceptable in very clean conditions	Within 4 days optimal 5–8 acceptable in very clean conditions	Within 6 months optimal Up to 12 months acceptable	Within 24 hours optimal
Thawed (previously frozen)	1–2 hours optimal	Up to 24 hours optimal	Never refreeze breast milk once thawed	Feed to infant within 24 hours
Leftover in bottle after feeding	Use within 2 hours of the start of your feed.			

Sources: Centers for Disease Control and Prevention (CDC), *Proper Storage and Preparation of Breast Milk*, November 27, 2023; Academy of Breastfeeding Medicine, *Human Milk Storage Information for Home Use for Full Term Infants*, 2017

BREAST MILK THAWING & SERVING GUIDELINES

Before serving, swirl the milk gently to incorporate any fat that has gathered on the side of the container. Milk that has been warmed (under running water or by your hand) will swirl more easily.

- Milk can be thawed in the refrigerator. It will take between 4 and 12 hours to defrost completely. Always thaw bags standing up and inside a clean container to avoid leaks.
- Milk can also be thawed under a stream of warm water, or by being placed in a bowl of warm water. If using a bowl, swirl gently to speed defrosting and to make sure that the fatty parts of the milk (which may separate when frozen) are incorporated. Make sure the temperature is no higher than 99 degrees. Never microwave breast milk!
- Fully defrosted milk can remain in the fridge for up to 24 hours, or at room temperature for two hours, before being used.
- If thawing milk still has ice crystals you can refreeze it, but any completely thawed milk that has been previously frozen should never be refrozen.
- Any breast milk that has been mixed with powdered or liquid formula should be used within an hour. Any milk left over should be disposed of.

Source: Centers for Disease Control and Prevention (CDC), Proper Storage and Preparation of Breast Milk, November 27, 2023

7. ROUTINE CLEANING

WARNINGS

Washable Parts

1. Do not sanitize Flange or Container parts using microwave or steam.
2. Do not use harsh chemicals to clean parts.
3. Do not store damp parts; fully air dry parts before storing.

Non-washable Parts

1. Do not sanitize Motor, Receiver, or Tubing using microwave or steam.
2. Never immerse Motor, Receiver, or Tubing in water.
3. Do not store Motor, Receiver, or Tubing in the refrigerator or freezer.

Safety


1. Wash hands with soap and water before handling.
2. Use a wash basin dedicated only to cleaning Pump parts and infant feeding items.

Note: For additional cleaning guidelines see the Center for Disease Control website, or your local health department guidance.

Washable Parts

Note: Avoid touching the inside of any parts that may come in contact with breast milk to prevent risk of contamination.


⚠️ Replace washable parts after signs of wear and tear.




Non-washable Parts

Motor, Receiver, and Tubing are NOT washable.
Wipe with disinfectant wipes (such as 70-90% ethyl or isopropyl alcohol).

⚠ Do not submerge or run under water.



Motor



Receiver



Tubing



NON-WASHABLE PARTS

CLEANING INSTRUCTIONS - AFTER EACH USE

Washable Parts

1. Disassemble parts.
2. Rinse milk-contacting parts immediately in warm water.
3. Soak with agitation for 5 minutes in hot, soapy water. Use mild dishwashing soap.
4. Wash each part for 1 minute in warm water.
5. Rinse well in warm water. If any residue is observed, repeat.
6. Air dry completely on a clean surface (e.g., clean paper towel or clean drying rack).
7. Store dry parts when not in use.

Non-washable Parts

1. Wipe with disinfectant wipes (such as 70-90% ethyl or isopropyl alcohol).

Storing the Willow Sync Breast Pump

Short Term

- Rinse milk-contacting parts in warm water.
- Soak with agitation for 5 minutes in hot, soapy water. Use mild dishwashing soap.
- Wash each part for 1 minute in warm water.

Long Term

- Thoroughly clean / dry washable parts.
- Prevent mold from forming by ensuring parts and storage bag are completely dry before storing.

8. CHARGING

- Check the battery status by briefly pressing the Power button.
- If the battery is low, charge with provided USB Charging Cable.
- When the battery is charging, the LED Indicator will pulse.
 - When the battery is fully charged, the LED Indicator light will stay solid Green.
- The battery can be fully charged in about 2 hours when using a 5V / 2A charger. Refer to page 25 for battery color indicators.

Important: Pump will not operate while plugged into electrical outlet. Only use included USB Charging Cable and chargers compatible with specifications identified on page 37 (Technical Information).

9. TROUBLESHOOTING

The Pump automatically turns off after 25 minutes. You may turn off the Pump at any time. Press the Power button to restart Pump with a new session.

If your Pump is not working, or the LED Indicator is blinking red:

- Make sure the battery is charged
- Make sure the Pump and Pump Cup are fully assembled
- Check parts for wear and tear; replace any damaged parts
- Reset by pressing and holding the Power button for at least 6 seconds

The Willow Sync™ Breast Pump is not a serviceable device. If you experience any issues or need help setting up or operating your Pump, please visit help.onewillow.com.

TECHNICAL INFORMATION

WARRANTY

Under our limited warranty, Willow * (“we,” “us,” “our”), warrants to the original consumer purchaser (the “purchaser”) to be free from defects in material and workmanship for the period of 12 months from the date of Pump purchase and the period of 90 days from purchase on all other product components under normal use and conditions (the “Warranty Period”). This limited warranty does not cover damage caused by accident, misuse, normal wear and tear, improper use, unauthorized maintenance or repair, or failure to follow written instructions for use. The Pump can not be submerged in water, washed or sanitized, placed in the refrigerator or freezer. The Flange can not be placed in microwave or microwave sterilizer bags. Doing so voids this limited warranty. Warranty coverage terminates if you sell or otherwise transfer this product to another person. This warranty gives you specific legal rights and you may also have other rights, which vary, by location.

WILLOW MAKES NO EXPRESS WARRANTY OF ANY KIND REGARDING THE PRODUCT OTHER THAN THOSE WARRANTIES EXPLICITLY SET FORTH HEREIN. ANY IMPLIED WARRANTY, INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY, NON-INFRINGEMENT, QUIET ENJOYMENT, OR ANY IMPLIED WARRANTY OF FITNESS FOR A PARTICULAR PURPOSE, TO THE EXTENT PERMITTED BY LAW, SHALL BE LIMITED IN DURATION TO A PERIOD OF 90 DAYS FROM THE DATE OF PURCHASE BY THE ORIGINAL PURCHASER. **Note: Willow Innovations, Inc.*

In the event the product is found by Willow to not meet the above limited warranty in the Warranty Period, as purchaser’s sole and exclusive remedy Willow will repair or, at the option of Willow, replace this product without charge for such replacement or parts or labor. NOTE: The purchaser shall bear all expense for returning this product to Willow. This limited warranty does not apply to any product that has been used commercially, subjected to misuse, abuse, or alteration.

Only Willow brand replacement parts are recommended for use with Willow breastfeeding products. Improper or incorrectly performed maintenance or repair, or use of non-Willow brand replacement parts may invalidate this limited warranty.

TO THE EXTENT PERMITTED BY LAW, Willow SHALL NOT BE LIABLE FOR ANY INDIRECT, SPECIAL, PUNITIVE, EXEMPLARY INCIDENTAL OR CONSEQUENTIAL DAMAGES INCLUDING, BUT NOT LIMITED TO, REPLACEMENT COSTS RESULTING FROM THE BREACH OF ANY WRITTEN OR IMPLIED WARRANTY.

If you wish to make a claim under this limited warranty, you must return the product to Willow prepaid, together with your dated bill of sale or other dated proof of purchase and a brief statement of the problem, along with a Returned Goods Authorization number. You may obtain the appropriate Authorization number and shipping address by calling the phone number. Items returned without an associated Returned Goods Authorization Number will be returned to sender.

To obtain the address to return goods or a Returned Goods Authorization Number contact the Customer Care Team.

We recommend recording your serial number here when you first open your Pump for easy reference when contacting Willow Customer Care Team. The serial number can be found on the bottom of the Pump.

Record the Serial Number(s) here for easy reference:

Motor											
-------	--	--	--	--	--	--	--	--	--	--	--

POWER

Use the USB Charging Cable included with your Pump.
USB charger is available at onewillow.com.

Recommended Charger Specifications:

Input: 100-240V AC, 50/60Hz 0.28A
Output: 5V DC 1.0A min. per port
Willow Charger Part Number:

REF WDCB76

Date of Manufacture:

Year Code																							
2024		2025		2026		2027		2028		2029		2030		2031									
I		J		K		L		M		N		O		P									
Month Code																							
Jan		Feb		Mar		Apr		May		Jun		Jul		Aug		Sep		Oct		Nov		Dec	
1		2		3		4		5		6		7		8		9		A		B		C	

SPECIFICATIONS

Pump dimensions: 80mm x 120mm x 135mm
Pump weight: 350g
Vacuum strength: 30-255 mmHg
Cycle speed: 20-100 cpm
Pump operating life: The operating life of this breast pump is defined to be approximately three 15-minute sessions per day, for 12 months. The operating life for the washable parts is 3 months.
Battery: Internal Lithium-Ion (not replaceable)

Date of Manufacture: The date of manufacture can be located within the serial number of the pump. The serial number is an 11 digit number formatted as xxxYMxxxxxx where Y is the year and M is the month. Use the tables to determine the year and month of manufacture.

OPERATING CONDITIONS

Operating Temperature	5 to 32 °C
Relative Humidity	15% – 90% non-condensing
Pressure	700 hPa to 1060 hPa

TRANSPORTATION AND STORAGE

Temperature	-20 to 60 °C
Relative Humidity	15% – 90% non-condensing
Pressure	700 hPa to 1060 hectopascal

DISPOSAL INFORMATION

The Willow Sync™ Breast Pump contains a rechargeable lithium ion battery. Use your local return and collection system for electrical waste and electronic equipment (including batteries). diaphragm, container assembly, and flange should be disposed of per your local waste management requirements.

ELECTROMAGNETIC COMPATIBILITY COMPLIANCE STATEMENT INFORMATION

The Willow Sync™ Breast Pump may be used at home, office, professional healthcare environment or other general areas.

Special consideration should be given to use:

- In High Electrostatic Discharge Environments: If the air is very dry, use the Willow Sync™ Breast Pump where floors are wood, concrete, or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%. If you notice shocks when you touch things, as sometimes happens in very dry winter months, be sure to touch a metal item such as a table before handling the Pump.

WARNINGS:

- Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.
- Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.
- Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the Willow Sync™ Breast Pump, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

The essential performance functions include vacuum not exceeding 295 mmHg and backflow protection to prevent milk leakage into the motor unit. The Willow Sync™ Breast Pump was tested for immunity to electromagnetic disturbances and passed using the criteria of maximum vacuum, backflow protection and operation of the breast pump.

Note: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

GUIDANCE AND MANUFACTURER'S DECLARATION - ELECTROMAGNETIC EMISSIONS

Emissions Test	Compliance	Electromagnetic Environment – Guidance
RF emissions CISPR 11	Group 1	The Willow Sync™ Breast Pump uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Group 1	The Willow Sync™ Breast Pump is intended for use in the electromagnetic environment specified below. The customer or user of the Willow Sync Breast Pump should ensure that it is used in such an environment.
RF emissions CISPR 11	Class B	The Willow Sync™ Breast Pump is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	Not applicable
Voltage fluctuations / flicker emissions IEC 61000-3-3	Complies	Not applicable

The Breast Pump has been tested and complies with the requirements for IEC 60601-1-2 4th Edition, Tables 4 through Table 9. See Guidance and Declaration of Electromagnetic Immunity in the following pages in this Instructions for Use.

The assembled Breast Pump is a Type BF Applied Part according to IEC 60601-1.

GUIDANCE AND MANUFACTURER'S DECLARATION - ELECTROMAGNETIC EMISSIONS

The Willow Sync™ Breast Pump is intended for use in the electromagnetic environment specified below. The customer or user of the Willow Sync Breast Pump should ensure that it is used in such an environment.

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment – Guidance
Electrostatic discharge (ESD) IEC 61000-4-2 Normal operation while on battery power and Battery Charge Mode	Contact: ±8 kV Air: ±2 kV, ±4 kV, ±8 kV, and ±15 kV	Contact: ±8 kV Air: ±2 kV, ±4 kV, ±8 kV, and ±15 kV	Floors should be wood, concrete, or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	±2 kV for power supply lines	±±2 kV for power supply lines	Mains power quality should be that of a typical commercial, home or hospital environment.
Surge IEC 61000-4-5	± 0.5 kV, ±1 kV line to line	± 0.5 kV, ±1 kV line to line	Mains power quality should be that of a typical commercial, home or hospital environment.

GUIDANCE AND MANUFACTURER'S DECLARATION – ELECTROMAGNETIC IMMUNITY

The Willow Sync™ Breast Pump is intended for use in the electromagnetic environment specified below. The customer or user of the Breast Pump should ensure that it is used in such an environment.

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment –Guidance
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	0 % UT; 0,5 cycle at 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0 % UT; 1 cycle at 0° 70 % UT; 25/30 cycles at 0° tested at 100 V 60 Hz and 240 VAC 50 Hz 0 % UT; 250/300 cycle tested at 230 V, 50 Hz	0 % UT; 0,5 cycle at 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0 % UT; 1 cycle at 0° 70 % UT; 25/30 cycles at 0° tested at 100 V 60 Hz and 240 VAC 50 Hz 0 % UT; 250/300 cycle tested at 230 V, 50 Hz	Mains power quality should be that of a typical commercial, home or hospital environment. If the user of the Willow Sync Breast Pump requires continued operation during power mains interruptions, it is recommended that the Willow Sync Breast Pump be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	30 A/m	30 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial, home or hospital environment.

NOTE: UT is the A.C. mains voltage prior to application of the test level.

GUIDANCE AND MANUFACTURER'S DECLARATION - ELECTROMAGNETIC EMISSIONS

The Willow Sync™ Breast Pump is intended for use in the electromagnetic environment specified below. The customer or user of the Willow Sync Breast Pump should ensure that it is used in such an environment.

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment –Guidance
Radiated RF EM Fields IEC 61000-4-3	10V/m 80MHz to 2.7 GHz	10V/m 80MHz to 2.7 GHz	Portable and mobile RF Communications equipment should be used no closer to any part of the Willow Sync Breast Pump, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance: $d = 1.2\sqrt{P}$ 80 MHz to 800 MHz 800 MHz to 2.7 GHz Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, - should be less than the compliance level in each frequency range. - Interference may occur in the vicinity of equipment marked with the following symbol: 
Conducted disturbances induced by RF fields IEC 61000-4-6	3 V; 150 kHz to 80 MHz 6 V; ISM and amateur radio bands between 0.15 MHz and 80 MHz	Not applicable	The Willow Sync™ Breast Pump is exempt from this test.

NOTE 1: At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.

ENCLOSURE PORT IMMUNITY TO RF WIRELESS COMMUNICATIONS EQUIPMENT

Test frequency (MHz)	Band (MHz)	Service	Modulation	Immunity test level (V/m)
385	380-390	TETRA 400	Pulse modulation 18 Hz	27
450	430-470	GMRS 460, FRS 460	FM ± 5 kHz deviation 1kHz sine	28
450	430-470	GMRS 460, FRS 460	FM ± 5 kHz deviation 1kHz sine	28
710	704-787	LTE Band 13, 17	Pulse modulation 217 Hz	9
745				
780				
810	800-960	GSM 800/900, TETRA 800, iDEN 820, CDMA 850, LTE Band 5	Pulse modulation 18 Hz	28
870				
930				
1720	1700-1990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1, 3, 4, 25; UMTS	Pulse modulation 217 Hz	28
1845				
1970				
2450	2400-2570	Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse modulation 217 Hz	28
5240	5100-5800	WLAN 802.11 a/n	Pulse modulation 217 Hz	9
5500				
5785				

SYMBOLS GLOSSARY (1 OF 2)

SYMBOL	DEFINITION	Standard from which the symbol was referenced, Reference/ clause/ registration number	SYMBOL	DEFINITION	Standard from which the symbol was referenced, Reference/ clause/ registration number
	Manufacturer	ISO 15223-1:2021, 5.1.1		Keep away from heat	ISO 15223-1:2021, 5.3.2
	Date of Manufacture	ISO 15223-1:2021, 5.1.3		Temperature Limit	ISO 15223-1:2021, 5.3.7
	Lot Number	ISO 7000, Reg. no 2492		Humidity Limitation	ISO 15223-1:2021, 5.3.8
	Catalog Number; also represents the model number	ISO 15223-1:2021, 5.1.6		Atmospheric Limitation	ISO 15223-1:2021, 5.3.9
	Serial Number	ISO 15223-1:2021, 5.1.7		Consult Instructions For Use	ISO 15223-1:2021, 5.4.3
	Do not use if package is damaged	ISO 7000, 2606		Direct Current	IEC 60417, 5031

SYMBOLS GLOSSARY (2 OF 2)

SYMBOL	DEFINITION	Standard from which the symbol was referenced, Reference/ clause/ registration number	SYMBOL	DEFINITION	Standard from which the symbol was referenced, Reference/ clause/ registration number
	Keep Dry	ISO 15223-1:2021, 5.3.4	IP22	Ingress Protection against entrance of foreign objects and against harmful effects due to the entrance of water.	IEC 60601-1, Table D.3 Symbol 2
	Class II Equipment – Device does not rely on earth plug or power cord.	IEC 60417, 5172		Do not throw in trash. Send to a recycling center	EN 50419, Marking of Electrical and Electronic Equipment in accordance with Article 11(2) of Directive 2002/96/EC (WEEE)
	Recycling	ISO 7000, 1135		Non-ionizing electromagnetic radiation present due to a RF transmitter.	IEC 60417, 5140
	Caution	ISO 15223-1:2021, 5.4.4		Product in compliance with UL's Listing and Follow-Up Service for Canada and the United States.	UL Classification mark
	Compliance with international requirements for protection from electric shock. (Type BF applied parts).	IEC 60601-1, Table D.1 Symbol			

FCC COMPLIANCE STATEMENT

CAUTION: Any changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate this equipment.

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions:

(1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

ISED STATEMENT

This device complies with Industry Canada license-exempt RSS standard(s).

Operation is subject to the following two conditions:

(1) this device may not cause interference, and (2) this device must accept any interference, including interference that may cause undesired operation of the device.



Manufactured for:

Willow Innovations, Inc
1975 W. El Camino Real #101
Mountain View, CA 94040
United States

Contact Information:

Customer Care
phone (USA): 1-888-WILLOW-3
web: help.onewillow.com

Designed in California, Made in Thailand

© 2025 Willow Innovations, Inc. Willow, Willow Go, Willow Sync, Willow 360, Flextube, the stylized mother baby logo and the stylized W logo are trademarks or registered trademarks of Willow in the United States and other countries.

Patent: onewillow.com/patents/



Medical - general medical equipment as to electrical shock, fire and mechanical hazards only in accordance with ANSI/AAMI ES60601-1(2005) + AMD 1(2012)
CAN/CSA-C22.2 No. 60601-1:14
E517883



CAREFULLY READ ALL INSTRUCTIONS PRIOR TO USE
Save instructions and sales receipt for future reference.



The Bluetooth® word mark and logos are registered trademarks owned by the Bluetooth SIG, Inc. and any use of such marks by Willow is under license. Other trademarks and trade names are those of their respective owners.

willowsync™