


medela




# Magic InBra™



**Before First Use** • Avant la première utilisation • Antes del primer uso •  
Före första användningen • Før første brug • Før første gangs bruk •  
Ennen ensimmäistä käyttöä






**Charge for 1 hour** • Chargez pendant  
1 heure • Cargar por 1 hora • Ladda i  
en timme • Oplad i 1 time • Lad i 1 time  
• Lataa 1 tunnin ajan



| <14 mm | 15 mm, 16 mm,<br>17 mm | 18 mm, 19 mm,<br>20 mm | 21 mm, 22 mm,<br>23 mm |
|--------|------------------------|------------------------|------------------------|
|        |                        |                        |                        |
|        |                        |                        |                        |

**Choose the right size** • Choisissez la bonne taille •  
Elige el tamaño correcto • Välj rätt storlek • Vælg den  
rigtige størrelse • Velg den  
riktige størrelsen • Valitse  
oikea koko





**Wash** • Lavez • Lava •  
Diska • Vask • Vask • Pesu

**Disinfect** • Désinfectez •  
Desinfecta • Desinficera •  
Desinficér • Desinfiser •  
Desinfiointi

**Dry** • Séchez • Seca •  
Torka • Tør • Tørk •  
Kuivaus

**Check assembly** • Vérifiez l'assemblage • Verifica el ensamblaje •  
Kontrollera montering • Kontrollér samlingen • Kontrollen montering •  
Tarkista kokoonpano



**After Pumping** • Après l'expression du lait • Después de la extracción de  
leche • Efter pumping • Efter pumpning • Etter pumping •  
Pumppaamisen jälkeen

